

Liquid Safety

SOLVENT MONITORINGSM-411 SERIES





SM

411

Level monitoring system for solvents, eluents and waste container

SM

411W

Including web-enabled remote monitoring via smartphone or tablet

BENEFITS

- All solvent levels of your HPLC system at a glance
- Avoids an empty solvent bottle as well as a spilling waste container
- Contactless level monitoring for highly pure eluents
- Warning in advance when a critical level is emerging
- Be informed of your solvent levels at any time at any place via the internet



PHONE: +41 61 413 80 42 EMAIL: info@hplconsult.ch



SOLVENT MONITORING

SM-411 SERIES

Technical Data

SM

411 / 411W



INPUTS: Up to 4 solvent/eluent bottles and 1 waste container
INDICATORS: 5 tricolor LEDs with traffic light principle: green, yellow & red

warning: Yellow signal when approaching a critical fill-level

• ALARM: Acoustic alarm and red signal after reaching the critical fill-level

output: Contact closures for solvent/eluent and waste levels

• configuration: Factory setting, customized settings possible

• SAFETY: Eliminate danger of waste overflow and running out of solvent/eluent

• communication: Fill-levels accessible via internet











PHONE: +41 61 413 80 42 **EMAIL**: info@hplconsult.ch